

# Signalpet x Provet Integration Setup guide

## Automatically push SignalPET reports into your Provet patient record

The SignalPET – Provet API integration allows SignalPET reports to be automatically pushed directly into the patient page in Provet. This helps keep all diagnostic insights centralized in one place, reduces manual uploads, and ensures your team can easily access reports within the patient's medical record.


## Choose which reports are sent to Provet

Veterinarians can control which types of SignalPET reports are included in the integration. You can select the report types that best fit your clinic's workflow, ensuring only relevant reports are automatically added to the Provet patient page.


# How to enable the integration in Provet

Enabling the SignalPET integration takes just a few clicks (screenshots from Provet below):

1. In Provet, go to Settings
2. Navigate to Integration
3. Click Add Application



4. Select the SignalPET application
5. Click Save



## What happens next?

Choose which reports are sent to Provet.

After the application is added, no additional setup is required from your side.

SignalPET will use the securely shared credentials to complete the integration setup for your account.

### Once enabled:

- SignalPET reports will be pushed directly into the relevant Provet patient page
- Your selected report types will appear automatically
- Your workflow remains uninterrupted, with no manual uploads needed

If you have any questions or would like help choosing which report types to enable, the SignalPET team is happy to assist.