

# SignalPET x Vetspire Integration Setup guide

## Automatically push SignalPET reports into your Vetspire patient record

The SignalPET – Vetspire API integration allows SignalPET reports to be automatically pushed directly into the patient page in Vetspire. This helps keep all diagnostic insights centralized in one place, reduces manual uploads, and ensures your team can easily access reports within the patient’s medical record.


## Choose which reports are sent to Vetspire

Veterinarians can control which types of SignalPET reports are included in the integration. You can select the report types that best fit your clinic’s workflow, ensuring only relevant reports are automatically added to the Vetspire patient page.


## How to enable the integration in Vetspire

Enabling the SignalPET integration takes just a few clicks (screenshots from Vetspire below):


1. In Vetspire, go to More
2. Navigate to Admin



3. Under Integrations section, click API Keys.




4. Click "Create New Key"



5. Fill the Title: "SignalPET API Key" Role: "API Key"

6. Click Create



## What happens next?

Choose which reports are sent to Vetspire.

After the API Key was created , add it to the integration setup page in SignalPET.

SignalPET will use the securely shared credentials to complete the integration setup for your account.

### Once enabled:

- SignalPET reports will be pushed directly into the relevant Vetspire patient page
- Your selected report types will appear automatically
- Your workflow remains uninterrupted, with no manual uploads needed

If you have any questions or would like help choosing which report types to enable, the SignalPET team is happy to assist.